



Biolmaging Systems | Ultraviolet Products | Laboratory Products | Light Sources

For information contact:  
Kathleen Buckman  
(909) 946-3197  
kbuckman@uvp.com

## Press Release

For Immediate Release

# New Digital Color Camera for Brilliant Color Images

March 2008

Upland, CA ... UVP, LLC has released the next generation digital color camera for the DigiDoc-It, ChromaDoc-It and MultiDoc-It series basic imaging systems. The Digi Camera 105 offers a sophisticated design with enhanced 12 megapixel resolution, the 6x optical zoom and high quality construction for researchers capturing fluorescent, colorimetric and TLC images. To quickly and easily acquire images, the Digi Camera 105 connects to your computer loaded with UVP's exclusive capture software. The software controls the camera image capture and allows you to see live preview of images as well as enhance and document images.



Each UVP system is uniquely designed to fit a variety of applications and budgets. The DigiDoc-It System is excellent for multiple user labs and features a 25-watt High Performance UV Transilluminator and lightweight hood. The MultiDoc-It Systems are versatile design with compact 8-watt Benchtop UV Transilluminator and sturdy hood for fluorescent imaging or no transilluminator and overhead UV lamps for TLC imaging. Add a UV to White Light Converter Plate to either system for colorimetric imaging of silver or coomassie blue stained images. ChromaDoc-It System is designed with overhead UV lamps specifically for TLC imaging applications.

For information, link to DigiDoc-It, MultiDoc-It, or ChromaDoc-It Systems.



**About UVP, LLC** - UVP is a global leader in life science imaging, supplying imaging systems for academic, biopharmaceutical and biotechnology applications. With core expertise as the leading developer and manufacturer of ultraviolet products since 1932, UVP continues to provide innovative products and specializes in fluorescence and luminescence-based imaging applications for the life sciences. The organization provides comprehensive service and support to customers and dealers worldwide. UVP's headquarters and manufacturing operations are located in Upland, California. European Operations, Ultra-Violet Products Ltd., is headquartered in Cambridge, England. For additional information, go to [www.uvp.com](http://www.uvp.com)

###